

LAPSIM® 2020 ANNUAL UPDATE

LAPAROSCOPY TRAINING JUST GOT BETTER

This is our biggest annual update yet! With a new-look user interface, improved feedback metrics and a new cloud portal so instructors can login online.

More Intuitive Interface

- An effective goal-orientated environment that helps residents learn and stay motivated.
- Easy to use and navigate, with easy-to-follow instructions and video guides.
- Clear overview of assigned courses.



Improved Metrics for Clear Objective Feedback

- Resident get a comprehensive data-driven training report.
- Assessments classified into easy to understand surgical terms, so residents know what they are good at and what they are not.
- All parameters are grouped according to the GOALS rating scale, the validated assessment tool that measures training outcomes.

New Cloud Portal

- Instructors can login anywhere, anytime.
- Manage admin quickly while on the go.



General Updates

- Maryland Jaw Vessel Sealer/Divider is now available with GYN/Salpingectomy.
- The energy instruments can now give residual heat in a more realistic manner - remember to be cautious and not cause burn damage after deactivating them!
- Plus a range of minor improvements and bug fixes.